

	<b>Information Sheet</b>	Rev. 0	P. 1 of 1
	<b>Instructions for the collection and storage of faeces for microbiological and parasitological examination</b>	P06_A0250_MAN001_FI001 Ref. P06_A0250_MAN001	

### CONTAINER

- Use the appropriate container with scoop
- If multiple investigations are required and the method is the same, a single container is sufficient.

### QUANTITY

- If the stool is solid, fill the scoop generously
- If the stool is liquid, fill the container up to the red line or fill at least 1/3 full.

### COPROCULTURE

**Store in a refrigerator and deliver within 24 hours if required for:**

- ⇒ GIARDIA SEARCH
- ⇒ CRYPTOSPORIDIUM
- ⇒ ANTIGENE H. PYLORI
- ⇒ C. DIFFICILE SEARCH

**Store in a refrigerator and deliver on the same day if requested for:**

- ⇒ VIRUS SEARCH (Rotavirus - Adenovirus - Norovirus)

### PARASITE SEARCH WITH MICROSCOPIC EXAMINATION

- It is advisable to collect three samples on alternate days
- Samples must be delivered in a single delivery if collected in bottles with preservative liquid

07/01/2024

Document owned by Territorial Social Health Company G. Pini - Orthopaedic Trauma Centre,  
for the exclusive use of recipients.

It may not be distributed in part or in whole to third parties without the written authorisation of the General Manager.  
Controlled distribution in electronic format.