

**Instructions for collection
and storage for faecal
occult blood testing**P06_A0250_MAN001_FI003
Ref. P06_A0250_MAN001

The Laboratory supplies 3 containers with screw caps containing preservative liquid

HOW TO COLLECT FAECAL SAMPLES:

- ⇒ There is no need to follow any special diets before the test;
- ⇒ Do not perform the test during your menstrual period or if you have haemorrhoids;
- ⇒ Avoid contaminating faeces with urine;
- ⇒ Collect stool samples on three different days (ideally on alternate days);
- ⇒ Unscrew and remove the GREEN cap from the container;
- ⇒ Drag the stick over the stool several times horizontally and vertically;
- ⇒ Insert the stick into the tube, screw it back on and shake;
- ⇒ Identify the samples using the labels provided with:
 - **Surname;**
 - **First Name;**
 - **Age;**
 - **Gender;**
 - **Collection date.**

Deliver the sample to the laboratory as quickly as possible.

If this is not possible, store in a refrigerator.

07/01/2024

Document owned by Territorial Social Health Company G. Pini - Orthopaedic Trauma Centre, for the exclusive use of the recipients.

It may not be distributed in part or in whole to third parties without the written authorisation of the General Manager.

Controlled distribution in electronic format.